


CABOS EXTENSORES

Oximetria - SpO²



EO-3181-0

EO-2181-1



EO	Conector Monitor SpO ²	Comp. cabo	Conector Sensor	Trava
	0 - Fio solto	06 - 60cm	0 - Fio solto	0 - Sem trava
	1 - DB7 Macho 	10 - 100cm	1 - DB9 Fêmea 	1 - Modelo BCI 
	2 - DB9 Macho 	18 - 180cm	////////////////////	2 - Modelo Nellcor 
	3 - Snap - 7 Pinos 	20 - 200cm	////////////////////	////////////////////
	////////////////////	30 - 300cm	////////////////////	////////////////////

Observações

1. Outros tipos de conectores para monitor sob consulta.
2. Outros Comprimentos de cabo sob consulta.
3. Pinos dos Contatos Torneados.

CATÉTER

Observações

Os cabos apresentam 2 metros como comprimento total padrão. Outras medidas, sob consulta.

Ponta Catéter




04



03

Cabo



00



04

CSW-03-04

↓
Ponta
Catéter

↘ Cabo

BBS

Eletromédicos LTDA.



Faça seu orçamento



Rua Florianópolis, 224 - São Paulo - SP -
CEP 03185-050



(+55) (11)2601-1051



WWW.BBSELETROMEDICOS.COM.BR



infotec@bbseletromedicos.com.br



EDIÇÃO: 08/05/2023

